



Recombinant Escherichia coli Glycerophosphoryl diester phosphodiesterase (ugpQ)

Product Code	CSB-EP320805ENV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10908
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSNWPYPRIV AHRGGGKLAP ENTLASIDVG AKYGHKMIEF DAKLSKDGEI FLLHDDNLER TSNGWGVAGE LNWQDLLRVD AGSWYSKMFK GEPLPLLSQV AERCREHGMM ANIEIKPTTG TGPLTGKMVA LAARELWAGM TPPLLSSFEI DALEAAQQAA PELPRGLLLD EWRDDWRELT ARLGCVSIHL NHKLLNKARV MQLKDAGLRI LVYTVNKPQR AAELLRWGVD CICTDAIDVI GPNFTAQ
Source	E.coli
Target Names	ugpQ
Protein Names	Recommended name: Glycerophosphoryl diester phosphodiesterase Short name= Glycerophosphodiester phosphodiesterase EC= 3.1.4.46
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.