



Recombinant Mouse mammary tumor virus Gag-Pro-Pol polyprotein (gag-pro-pol), partial

| | |
|--------------------------|---|
| Product Code | CSB-EP320849MIU |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P11283 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mouse mammary tumor virus (strain C3H) (MMTV) |
| Purity | >85% (SDS-PAGE) |
| Source | E.coli |
| Target Names | gag-pro-pol |
| Protein Names | Recommended name: Gag-Pro-Pol polyprotein Cleaved into the following 10 chains: 1. Matrix protein p10 2. Phosphorylated protein pp21 3. Protein p3 4. Protein p8 5. Protein n 6. Capsid protein p27 7. Nucleocapsid protein-dUT |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |