



Recombinant *Saccharomyces cerevisiae* Nuclear fusion protein BIK1 (BIK1)

Product Code	CSB-MP320897SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11709
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDRYQRKIGC FIQIPNLGRG QLKYVGPVDT KAGMFAGVDL LANIGKNDGS FMGKKYFQTE YPQSGLFIQL QKVASLIEKA SISQTSRRTT MEPLSIPKNR SIVRLTNQFS PMDDPKSPTP MRSFRITSRH SGNQQSMDQE ASDHHQQQEF GYDNREDRME VDSILSSDRK ANHNTTSDWK PDNGHMNDLN SSEVTIELRE AQLTIEKLQR KQLHYKRLLD DQRMVLEEVQ PTFDRYEATI QEREKEIDHL KQLELELERRQ QAKQKQFFDA ENEQLLAVVS QLHEEIKENE ERNLSHNQPT GANEDVELLK KQLEQLRNIE DQFELHKTW AKEREQLKMH NDSLSKEYQN LSKELFLTKP QDSSSEEVAS LTKKLEEANE KIKQLEQAQA QTAVESLPIF DPPAPVDTTA GRQQWCEHCD TMGHNTAECF HHNPDNQQFF
Source	Mammalian cell
Target Names	BIK1
Protein Names	Recommended name: Nuclear fusion protein BIK1
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.