



Recombinant Rat Numb-like protein (NumbI)

Product Code	CSB-MP016187RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1L1I3
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MSRSAAASGG PRRPDQHLPP APCGASGPPE TFRTESDGAG TMNKLRQSLR RRPAYVPEA SRPHQWQADE DAVRKGTCSE PVRYLGHVEV EESRGMHVCE DAVKCLKAMG RKS VKSVLWV SADGLRVVDD KTKDLLVDQT IEKVSFCAPD RNLDKAFSYI CRDGTTRRWI CHCFLALKDS GERLSHAVGC AFAACLERKQ RREKECGVTA AFDA SRTSFA REGSFRLSGG GRPAEREAGD KKKAEAAAAP AVAPGPAQPG HVSPTPATT PGEKGEAGTP VAAGTTAAAI PRRHAPLEQL VRQGSFRGFP ALSQKNSPFK RQLSLRLNEL PSTLQRRTDF QVKGTVPEME PPGTGSDSDGI SALCTQISS FASAGAPVPG PPSATTGTSA WGEPSVPAAT AFQPGHKRTP SEAERWLEEV SQVAKAQQQQ QQQQQQQQQQ QQQQQQQQQQ QATSVPPMPT MAPTLQPFST PVGPFDTAAA QVAVFLPPTH MQPPFVPAYP GLGYPPMPRV PVVGITPSQM VANAFCSAAQ LQPQPATLLG KAGAFPPTA PSAPGGQARP RPNGAPWPPE PAPAPAPELD PFEAQWAALE GKPAVEKPSN PFSGDLQKTF EIEL
Source	Mammalian cell
Target Names	NumbI
Protein Names	Recommended name: Numb-like protein
Expression Region	1-614
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.