



Recombinant Human RNA-binding protein MEX3A (MEX3A)

| | |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP013740HU-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | A1L020 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MPSLVVSGIM ERNGGFGELG CFGGSAKDRG LLEDERALQL ALDQLCLLGL GEPAPTAGE DGGGGGGGAP AQPAAPPQPA PPPPPAAPPA APTAAPAAQT PQPPTAPKGA SDAKLCALYK EAELRLKGSS NTTECVPVPT SEHVAEIVGR QGCKIKALRA KTNTYIKTPV RGEFPVFMVT GRREDVATAR REIISAAEHF SMIRASRNKS GAAFGVAPAL PGQVTIRVRV PYRVVGLLVG PKGATIKRIQ QQTNTYIITP SRDRDPVFEI TGAPGNVERA REEIETHIAV RTGKILEYNN ENDFLAGSPD AAIDSRYSDA WRVHQPGCKP LSTFRQNSLG CIGECGVDSG FEAPRLGEQG GDFGYGGYLF PGYGVGKQDV YYGVAETSPP LWAGQENATP TSVLFSSASS SSSSSAKARA GPPGAHRSPA TSAGPELAGL PRRPPGEPLQ GFSKLGGGGL RSPGGGRDCM VCFESEVTAA LVPCGHNLFM MECAVRICER TDPECPVCHI TATQAIRIFS |
| Source | E.coli |
| Target Names | MEX3A |
| Protein Names | Recommended name: RNA-binding protein MEX3A Alternative name(s): RING finger and KH domain-containing protein 4 |
| Expression Region | 1-520 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |