



Recombinant Drosophila melanogaster Enhancer of split m4 protein (m4)

Product Code	CSB-MP321060DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P13095
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MCQNKINTNN TMTIKSNKKL SYSVKLLQK IFKQQQRVEE EQNLKNALKA NSLESLESME NSRNADLESA SICASLESCE NEANERLSQS CEIEDYDFEQ LPTVPVHFVR TAHGTFFWTA VSDLPADNDL VEPLYCSTSN AIAIPQDRWV QA
Source	Mammalian cell
Target Names	m4
Protein Names	Recommended name: Enhancer of split m4 protein Short name= E(spl)m4
Expression Region	1-152
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.