



Recombinant Human Protein FAM22F (FAM22F)

Product Code	CSB-MP008184HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1L443
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASNGAYPVL GPGVTVNPGT SLSVFTALPF ATPAPGPAHR PPLVTAVVPP AGPLVLSAFP STPLVAGQDG RGPSEGAGASN VFVQMRTEVG PVKPPQAQTL ILTQAPLVWQ APGTLCGGVM CPPPLLLAAA PGVPVTSQV VGGTQACEGG WSHGLPLPPP PPAAQVAPIV SPGNARPWPQ GAHGEGLAP SQAKARPDDS CKPKSVYENF RLWQHYKPLA RRHLPQSPDT EALSCFLIPV LRSLARRKPT MTLEELGWQA MREWQHTSNF DRMIFYEMAE KFLEFEAEAE MQIQKSQWMK GPQSLPPPAP PRLEPRGPPA PEVVKQPVYL PSKDGPKAPT ACLPPPRPQR PAETKAHLPP PRPQRPAETN AHLPPPRPQR PAETKVPEEI PPEVVQEYVD IMEELLGSHP GDTGEPEGQR EKGKVEQPQE EDGITSDPGL LSYIDKLCSSQ EDFVTKEAV IHPRFLEELL SPDPQMDFLA LSQELEQEEG LTLAQLVEKR LLSLKEKGC RAAPRHGTAR LDSSPSEFAA GQEAAREVPD PQQRVSVETS PPQTAAQDPQ GQGRVRTGMA RSEDPAVLLG CQDSPRLKAV RPTSPQDHR PTCPLGLTKD ALGLPGESPV KESHGLAKGS SEETELPGMV YVVGSHHRLR PWRLSQSPVP SSGLLSPGGR GPQGALQSPS AQKRGLSPSP SPASKSKKRP LFGSPSPA EK TPHPGPGLRV SGEQSLAWGL GGPSQSQRK GDPLASRRKK KRHCSQ
Source	Mammalian cell
Target Names	NUTM2F
Protein Names	Recommended name: Protein FAM22F
Expression Region	1-756
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.