



Recombinant Lolium perenne Pollen allergen Lol p 1

Product Code	CSB-MP321329LQK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P14946
Product Type	Recombinant Protein
Immunogen Species	Lolium perenne (Perennial ryegrass)
Purity	>85% (SDS-PAGE)
Sequence	IAKVPPG PNITAEYGDK WLDAKSTWYG KPTGAGPKDN GGACGYKNVD KAPFNGMTGC GNTPIFKDGR GCGSCFEIKC TKPESCSGEA VTVTITDDNE EPIAPYHFDL SGHAFGSMK KGEEQNVRSA GELELQFRRV KCKYPDDTKP TFHVEKASNP NYLAILVKYV DGDGDVVAVD IKEKGKDKWI ELKESWGAVW RIDTPDKLTG PFTVRYTTEG GTKSEFEDVI PEGWKADTSY SAK
Source	Mammalian cell
Protein Names	Recommended name: Pollen allergen Lol p 1 Alternative name(s): Allergen Lol p I Allergen R7 Allergen= Lol p 1
Expression Region	24-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.