



Recombinant Escherichia coli Fe (3+) dicitrate transport ATP-binding protein FecE

Product Code	CSB-MP321341ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15031
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MTLRTENLTV SYGTDKVLND VLSLPTGKI TALIGPNGCG KSTLLNCFSR LLMPQSGTVF LGDNPINMLS SRQLARRLSL LPQHHLTPEG ITVQELVSYG RNPWLSLWGR LSAEDNARVN VAMNQTRINH LAVRRLTELS GGQRQRAFLA MVLAQNTPVV LLDEPTYLD INHQVDLMRL MGELRTQGKT VVAVLHDLNQ ASRYCDQLVV MANGHVMAQG TPEEVMTPGL LRTVFSVEAE IHPEPVSGRP MCLMR
Source	Mammalian cell
Target Names	fecE
Protein Names	Recommended name: Fe(3+) dicitrate transport ATP-binding protein FecE Alternative name(s): Iron(III) dicitrate transport ATP-binding protein FecE
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.