



Recombinant Zea mays Malate dehydrogenase [NADP], chloroplastic

Product Code	CSB-BP321391ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15719
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	SVDATKQAQD GVATAVATEA PASRKECFGV FCTTYDLKAE DKTKSWRKLV NVAVSGAAGM ISNHLLFKLA SGEVFGQDQP IALKLLGSER SFQALEGVAM ELEDSELYPLL REVSIGIDPY VVFQDVDWAL LIGAKPRGPG MERAALLDIN GQIFADQGKA LNAVASRNDE VLVVGNPCNT NALICLKNAP NIPAKNFHAL TRLDENRAKC QLALKAGVFY DKVSNVTIWG NHSTTQVPDF LNAKIDGRP V KEVIKDTKWL EEEFTLTVQK RGGVLIQKWG RSSAASTAVS IVDAIRSLVT PTPEGDFWST GUYTTGNPYG IAEDIVFSMP CRSKGDGDYE LASDVLMDDF LWERIKKSEA ELLAEKKCVA HLTGEGNAFC DLPEDTMLPG EV
Source	Baculovirus
Protein Names	Recommended name: Malate dehydrogenase [NADP], chloroplastic EC=1.1.1.82 Alternative name(s): NADP-MDH
Expression Region	41-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.