




Recombinant Streptomyces hygroscopicus Phosphinothricin N-acetyltransferase (bar)

Product Code	CSB-YP321441SOO
Relevance	This enzyme is an effector of phosphinothricin tripeptide (PTT or bialaphos) resistance. Inactivates PTT by transfer of an acetyl group.
Abbreviation	Recombinant Streptomyces hygroscopicus bar protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16426
Alias	Phosphinothricin-resistance protein
Product Type	Recombinant Protein
Immunogen Species	Streptomyces hygroscopicus
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MSPERRPADIRRATEADM PAVCTIVNHYIETSTVNFRTPEPQEPQEWTDLLVRL RERYPWLVAEVDGEVAGIAYAGPWKARNAYDWTAESTVYVSPRHQRTGLGS TLYTHLLKSLEAQGFKSVVAVIGLPNDPSVRMHEALGYAPRGMLRAAGFKHGN WHDVGFWQLDFSLPVPVPRPVLVPTETI
Research Area	Others
Source	Yeast
Target Names	bar
Protein Names	Recommended name: Phosphinothricin N-acetyltransferase Short name= PPT N-acetyltransferase EC= 2.3.1.183 Alternative name(s): Phosphinothricin- resistance protein
Expression Region	1-183aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	22.6kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.