



Recombinant Mouse TANK-binding kinase 1-binding protein 1 (Tbkbp1)

Product Code	CSB-EP023233MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2A9T0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MESMFEDDIS ILTQEALGPS EVWLDGPGDP SLGGDMCSAS HFALITAYGD IKERLGGLER ENATLRRRLK VYEIKYPLIT DFGEEHGFPL YELKDGSLL VEKVSLQQRL NQFQHELQKS KEQEEQLGEM IQAYEKLCVE KSDLETELGE MRALVETHLR QICGLEKQLQ QQQGLRDAAF SSLSPPAVPA SACPDLDLHY LALRGGPALG HAGWPGPSTV SVSELERRL EEALEAAQGE ARGAQLREEQ LQAECERLQG ELKQLQETRA QDLASNQSEC DMAWVKRVGD DQVNLALAYT ELTEELGRLR ELSSLQGRIL RTLLQEARN AGQRHSPLSQ RHSPAPACPS PSPPARPPPC APCQSPAAQR RSPVPPCPSP QRRRSPASPS CPSPVPQRRS PVPPSCQSPS PRRRSPVPPS CPAPQPRPPP PPGERTLAER VYAKPPSHA KAGFQGRRSY SELAEGAAYA GASPAWLQAE AATLPKPRAY GGELYGRPLS PRRAFEGIRL RFEKQPSEEE EWAMPASPPS PEAGTIRCAS FCAGFPIPES PAATAYAHAE HAQSWPSINL LMETVGS DIR SCPLCQLGFP VGYPDDALIK HIDSHLENSK I
Source	E.coli
Target Names	Tbkbp1
Protein Names	Recommended name: TANK-binding kinase 1-binding protein 1 Short name= TBK1-binding protein 1 Alternative name(s): Similar to NAP1 TBK1 adapter
Expression Region	1-611
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.