



Recombinant *Saccharomyces cerevisiae* 60S ribosomal protein L8-A (RPL8A)

Product Code	CSB-BP321500SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17076
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	APGKKVAPA PFGAKSTKSN KTRNPLTHST PKNFGIGQAV QPKRNLSRYV KWPEYVRVQR QKKILSIRLK VPPTIAQFQY TLDNRNTAAET FKLFNKYRPE TAAEKKERLT KEAAVAEGK SKQDASPKPY AVKYGLNHVV ALIENKKAKL VLIANDVDPI ELVVFLPALC KKMGGVPYAIV KGKARLGTLV NQKTSAVAAL TEVRAEDEAA LAKLVSTIDA NFADKYDEVK KHWGGGILGN KAQAKMDKRA KNSDSA
Source	Baculovirus
Target Names	RPL8A
Protein Names	Recommended name: 60S ribosomal protein L8-A Alternative name(s): L4 L4-2 L7a-1 Maintenance of killer protein 7 RP6 YL5
Expression Region	2-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.