



Recombinant Avian erythroblastosis virus Transforming protein Ski (V-SKI)

Product Code	CSB-EP321561ARH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17863
Product Type	Recombinant Protein
Immunogen Species	Avian erythroblastosis virus (strain Sloan-Kettering)
Purity	≥85% (SDS-PAGE)
Sequence	FHLSSMSSLG GPAAFSARWA QEMYKKDNGK DPAEPVLHLP PIQPPPVMPG PFFMPSDRST ERCETILEGE TISCFVVGGE KRLCLPQILN SVLRDFSLQQ INSVCDELHI YCSRCTADQL EILKVMGILP FSAPSCGLIT KTD AERLCNA LLYGGTYPPH CKKEFSSTIE LELTEKSFKV YHECFGKCKG LLVPELYSNP SAACIQCLDC RLMYPPHKFV VSHHSLENR TCHWGFDSAN WRSYILLSQD YTGKEEKARL GQLLDEMKEK FDYNNKYKRK APRNRESRPRV QLRRNKMFKT MLWDPAGGSA VLQRQPDGNE VPSDPPASKK TKIDDSASQS PASTEKEKQS SRLRSLSSSS NKSIGCVHPR QRLSAFRPWS PAVSANEEKEL STHLPALIRD SSFYSYKSFE NAVAPNVALA PPAQQKVVS NPPCATVV
Source	E.coli
Target Names	V-SKI
Protein Names	Recommended name: Transforming protein Ski
Expression Region	1-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.