



Recombinant Escherichia coli Uncharacterized protein yfaA (yfaA)

Product Code	CSB-EP321574ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17994
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSGKKAKGW RFYGLVGFGA IALLSAGVWA LQYAGSGPEK TLSPLVVHNN LQIDLNEPDL FLDSDSLSQL PKDLLTIPFL HDVLSDFVF YYQNHADRLG IEGSIRRIVY EHDLTLDKDL FSSLLDQPAQ AALWHDKQGH LSHYMVLIQR SGLSKLLEPL LFAATSDSQL SKTEISSIKI NSETVPVYQL RYNGNNALMF ATYQDKMLVF SSTDMLFKDD QQDTEATAIA GDLLSGKKRW QASFGLEERT AEKTPVRQRI VVSARWLGFG YQRLMPSFAG VRFEMGNDGW HSFVALNDES ASVDASFDFT PVWNSMPAGA SFCVAVPYSH GIAEEMLSHI SQENDKLNGA LDGAAGLCWY EDSKLQTPLF VGQFDGTAEQ AQLPGKLFTQ NIGAHESKAP EGVLPVSQTQ QGEAQIWRRE VSSRYGQYPK AQAAQPDQLM SDYFFRVSLA MQNKTLFSL DDTLVNNAALQ TLNKTRPAMV DVIPTDGIVP LYINPQGIK LLRNETLTSL PKNLEPVFYN AAQTLLMPKL DALSQQPRYV MKLAQMEPGA AWQWLPITWQ PL
Source	E.coli
Target Names	yfaA
Protein Names	Recommended name: Uncharacterized protein yfaA
Expression Region	1-562
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.