



Recombinant Replication initiation protein (repZ)

Product Code	CSB-YP321577ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P18023
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	≥85% (SDS-PAGE)
Sequence	MAGLKNTPIN AVHWSQLAPE EQIRFWEDYE AGRATTFLVE PERKRTKRRR GEHSTKPKCE NPSWYRPERY KALKGQLGHA YNRLVKKDPV TGEQSLRMRM SRHPFYVQKR TFVGRKYAFR PEKQRLLDAL WPVLVSFSDA GTHTVGMSVT RLAEIISPDK SEGHVIPELE VTVSRLSRL AEQVRFGVLG VSEETMWDRE HRQRLPRYVW ITPAGWQMLG VDMVKLHEQQ QKRLRESEIR QQLIREGVLR EDEDISVHAA RKRWYLQRSQ DALKKRREKA AASKRANRLK KLPVDQQIYE MAEYLRKRLP PDEAYFCSDD HLRRLAIREL RQLELTLAAP PPH
Source	Yeast
Target Names	repZ
Protein Names	Recommended name: Replication initiation protein Alternative name(s): Replication-associated protein
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.