



Recombinant *Saccharomyces cerevisiae* DNA-directed RNA polymerases I, II, and III subunit RPABC3 (RPB8)

Product Code	CSB-BP321805SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20436
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	SNTLFDDIF QVSEVDPGRY NKVCRIEAAAS TTQDQCKLTL DINVELFPVA AQDSLTVTIA SSLNLEDTPA NDSSATRSWR PPQAGDRSLA DDYDYVMYGT AYKFEEVSKD LIAVYYSFGG LLMRLEGNYR NLNNLKQENA YLLIRR
Source	Baculovirus
Target Names	RPB8
Protein Names	Recommended name: DNA-directed RNA polymerases I, II, and III subunit RPABC3 Short name= RNA polymerases I, II, and III subunit ABC3 Alternative name(s): ABC14.4 ABC14.5 DNA-directed RNA polymerases I, II, and III 14.5 kDa polypepti
Expression Region	2-146
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.