



# Recombinant Emericella nidulans Acyl-coenzyme A:6-aminopenicillanic-acid-acyltransferase 40 kDa form (penDE)

<b>Product Code</b>	CSB-EP321895EKF-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P21133
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Emericella nidulans (strain FGSC A4 / ATCC 38163 / CBS 112.46 / NRRL 194 / M139) (Aspergillus nidulans)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLHVTCQGTP SEIGYHHGSA AKGEIAK AID FATGLIHGKT KKTQAELEQL LRELEQVMKQ RWPRYYEEIC GIAKGAEREV SEIVMLNTRT EFAYGLVEAR DGCTTVYCKT PNGALQGQNW DFFTATKENL IQLTICQPGL PTIKMITEAG IIGKVGFN SA GVAVNYNALH LHGLRPTGLP SHLALRMALE STSPSEAYEK IVSQGGMAAS AFIMVGN AHE AYGLEFSPIS LCKQVADTNG RIVHTNHCLL NHGPSAQELN PLPDSWSRHG RMEHLLSGFD GTKEAFAKLW EDEDNYPLSI CRAYKEGKSR GSTLFNIVFD HVGRKATVRL GRPNNPDETF VMTFSNLDTK SAIQANI
<b>Source</b>	E.coli
<b>Target Names</b>	penDE
<b>Protein Names</b>	Recommended name: Acyl-coenzyme A:6-aminopenicillanic-acid-acyltransferase 40 kDa form EC= 2.3.1.164 Alternative name(s): Isopenicillin-N N-acyltransferase Cleaved into the following 2 chains: 1. Acyl-coenzyme A:6-aminopenicillanic-acid
<b>Expression Region</b>	1-357
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.