



Recombinant Retron EC67 protein

Product Code	CSB-MP321935ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21325
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MTKTSKLDAL RAATSREDLA KILDIKLVFL TNVLYRIGSD NQYTQFTIPK KGKGVRTISA PTDRLKDIQR RICDLLSDCR DEIFAIRKIS NNYSFGFERG KSIILNAYKH RGKQIILNID LKDDFFESFNF GRVRYGFLSN QDFLLNPVVA TTLAKAACYN GTLPQGSPCS PIISNLICNI MDMRLAKLAK KYGCTYSRYA DDITISTNKN TFPLEMATVQ PEGVVLGKVL VKEIENSGFE INDSKTRLTY KTSRQEV TGL TVNRIVNIDR CYYKKTRALA HALYRTGEYK VPDENGVLVS GGLDKLEGMF GFIDQVDKFN NIKKKLNKQP DRYVLTNATL HGFKLKLNAR EKAYSKFIYY KFFHGNTCPT IITEGKT DRI YLKAALHSLE TSYPELFREK TDSKKKEINL NIFKSNEKTK YFLDLGGTA DLKKFVERYK NNYASYGGSV PKQPVIMVLD NDTGPSDLLN FLRNKVKSCP DDVTEMRKMK YIHVFYNYLI VLTPLSPSGE QTSMEDLFPK DILDIKIDGK KFNKNNDGDS KTEYGKHIFS MRVVRDKKRK IDFKAFCCIF DAIKDIKEHY KLMLNS
Source	Mammalian cell
Protein Names	Recommended name: Retron EC67 protein Alternative name(s): ORF4-EC67 RT Including the following 2 domains: RNA-directed DNA polymerase EC= 2.7.7.49 Alternative name(s): Reverse transcriptase Ribonuclease H Short name=
Expression Region	1-586
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.