



Recombinant Mouse tRNA pseudouridine synthase-like 1 (Pus11)

Product Code	CSB-MP019091MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2ADA5
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MGSCGAVGSV RARYLVFLQY LGTDFNGVAA VRGNPRAVGV LNFLEEAAKR LNSVDPVRFT ISSRTDAGVH ALSNA AHLDI QRRPGQSPFS PEVVAKALNT HLKHPAIRVL KAFRVPNDFH ARHAATSRTY QYRLATGCSW PNQLPVFEQN VCWALQTEYL DMAAMQEAAQ HLLGTHDFSA FQSAGSPVTN TVRTLRRVSV SPGPASPFVL PEGSRRLQFW TLEFESQSFL YRQVRRMTAV LVAVGLGILA PTQVKVILES QDPLGKYQAR VAPARGLFLK SVLYDNFGPT S
Source	Mammalian cell
Target Names	Pus11
Protein Names	Recommended name: tRNA pseudouridine synthase-like 1 EC= 5.4.99.- Alternative name(s): tRNA pseudouridylate synthase-like 1 tRNA-uridine isomerase-like 1
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.