



Recombinant Bacillus subtilis 30S ribosomal protein S2 (rpsB)

Product Code	CSB-MP321964BRJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21464
Product Type	Recombinant Protein
Immunogen Species	Bacillus subtilis (strain 168)
Purity	>85% (SDS-PAGE)
Sequence	SVISMKQLL EAGVHFGHQT RRWNPCKMKRY IFTERNGIYI IDLQKTVKKV EEAYNFTKNL AAEGGKILFV GTKKQAQDSV KEEAQRSGMY YVNQRWLGGT LTFETIQKR IKRLKDIEKM QENGTDFVLP KKEVVQLKKE LERLEKFLGG IKDMKDLFDA LFIIDPRKER IAVAEARKLN IPIIGIVDTN CDPDEIDVVI PANDDAIRAV KLLTSKMADA ILEAKQGEEE AEVAEETAPE TETTA
Source	Mammalian cell
Target Names	rpsB
Protein Names	Recommended name: 30S ribosomal protein S2 Alternative name(s): BS1 Vegetative protein 209 Short name= VEG209
Expression Region	2-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.