



Recombinant *Synechococcus elongatus* 6-phosphogluconate dehydrogenase, decarboxylating (gnd)

Product Code	CSB-BP321986FPY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21577
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MALQQFGLIG LAVMGENLAL NIERNGFSLT VYNRTAEKTE AFMADRAQ GK NIVPAYSLED FVASLERPRR ILMVKAGGP VDAVVEQLKP LLDPGDLIID GGNSLFTDTE RRVKDLEALG LGFMGMGVSG GEEGALNGPS LMPGGTQAAY EAVEPIVRSI AAQVDDGPCV TYIGPGGSGH YVKMVHNGIE YGDMQLIAEA YDLLKSVAGL NASELHDVFA AWNKTPELDS FLIEITADIF TKVDDLGTGQ PLVELILDAA GQKGTGRWTV ETALEIGVAI PTIAAVNAR ILSSIKAERQ AASEILSGPI TEPFSGDRQA FIDSVRDALY CSKICSYAQQ MALLAKASQV YNYGLNLGEL ARIWKGGCII RAGFLNKIKQ AYDADPTLAN LLLAPEFRQT ILDRQLAWRR VIAIAAERGI PVPAFSASLD YFDSYRRDRL PQNLTQAQRD YFGAHTYERT DRSGSFHAQW F
Source	Baculovirus
Target Names	gnd
Protein Names	Recommended name: 6-phosphogluconate dehydrogenase, decarboxylating EC= 1.1.1.44
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.