



Recombinant *Pantoea ananas* Lycopene cyclase (crtY)

Product Code	CSB-MP322008ERA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21687
Product Type	Recombinant Protein
Immunogen Species	<i>Pantoea ananas</i> (<i>Erwinia uredovora</i>)
Purity	>85% (SDS-PAGE)
Sequence	MQPHYDLILV GAGLANGLIA LRLQQQPDM RILLIDAAPQ AGGNHTWSFH HDDLTESQHR WIAPLVVHHW PDYQVRFPTTR RRKLNSGYFC ITSQRFAEVL QRQFGPHLWM DTAVAEVNAE SVRLKKGQVI GARAVIDGRG YAANSALSVG FQAFIQEWR LSHPHGLSSP IIMDATVDQQ NGYRFVYSLP LSPTRLLIED THYIDNATLD PECARQNICD YAAQQGWQLQ TLLREEQGAL PITLSGNADA FWQQRPLACS GLRAGLFHPT TGYSLPLAVA VADRLSALDV FTSASIHAI THFARERWQQ QGFFRMLNRM LFLAGPADSR WRVMQRFYGL PEDLIARFYA GKLTLDRLR ILSGKPPVPV LAALQAIMTT HR
Source	Mammalian cell
Target Names	crtY
Protein Names	Recommended name: Lycopene cyclase
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.