



Recombinant *Saccharomyces cerevisiae* Ubiquitin-conjugating enzyme E2 1 (UBC1)

Product Code	CSB-EP322016SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21734
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSRAKRIMKE IQAVKDDPAA HITLEFVSES DIHHLKGTFL GPPGTPYEGG KFVVDIEVPM EYPFKPPKMV FDTKVYHPNI SSVTGAICLD ILKNAWSPVI TLKSALISLQ ALLQSPEPND PQDAEVAQHY LRDRESFNKT AALWTRLYAS ETSNGQKGNV EESDLYGIDH DLIDEFESQG FEKDKIVEVL RRLGVKSLDP NDNNTANRIL EELLK
Source	<i>E.coli</i>
Target Names	UBC1
Protein Names	Recommended name: Ubiquitin-conjugating enzyme E2 1 EC= 6.3.2.19 Alternative name(s): Ubiquitin carrier protein Ubiquitin-conjugating enzyme E2-24 kDa Ubiquitin-protein ligase
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.