



Recombinant *Phytophthora infestans* Actin-2 (ACTB)

Product Code	CSB-YP322091PJO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22132
Product Type	Recombinant Protein
Immunogen Species	<i>Phytophthora infestans</i> (Potato late blight fungus) (<i>Botrytis infestans</i>)
Purity	>85% (SDS-PAGE)
Sequence	MDDDIQAVVI DNGSGMCKAG FAGDDAPRAV FPSIVGMPKH LGIMVGMNQK DAYIGDEAQA KRGVLTLYP IEHGIVTNWD DMEKIWSHTF YNELRVAPEE HPVLLTEAPL NPKANRERMT QIMFETFNVP AMYVNIQAVL SLYASGRRTTG CVLDSGDGVV HTVPIYEGYA LPHAIVRDL AGRDLTDYMM KILTERGYSF TTTAEREIVR DIKEKLTIVA MDFDEEMEKS TRSSALDKTY ELPDGNVIVI GNERFRTPEV LFNPSMIGRE CSGVHECAFQ TIMKCDVDIR RDLYSNVVLS GGSTMFPGIG ERMTKEVIKL APTAMKVKII TPPERKYSVW IGGSILASLA TFQHMWISKT EYDESGPSIV HRKCF
Source	Yeast
Target Names	ACTB
Protein Names	Recommended name: Actin-2
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.