



Recombinant *Drosophila melanogaster* Protein embryonic gonad (eg)

Product Code	CSB-YP322205DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15370
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MNQLCKVCGE PAAGFHFGAF TCEGCKSFFG RTYNNIAAIA GCKHNGDCVI NKKNR TACKA CRLRKCLLVG MSKSGSRYGR RSNWFKIHCL LQEQQTTSGL GGSSVSGSGS GGGVSSASLE QLARLQQASN QARQTYQDKT NPCIKSATAT TSPRIEGAAV GTGIGGGASP SFLQAAKLHH QRQLKLSRL SNTPSDSGAS SAGDPNEDGV TSVLGGQIAT PSSTNATSLP KLDLRHPNFP ATSEPDADMQ RQRHQELLEI FRSHSEPLYS SFAPFSLP VLLAAGVPQL PIFKDQFKAE LLFPTTSSPE LEEPIDLSFR SRADPASPMA HNSNSPSLSE PAAASHCLGE STNFVRKSTP LDLTLVRSQT LTG
Source	Yeast
Target Names	eg
Protein Names	Recommended name: Protein embryonic gonad Alternative name(s): Nuclear receptor subfamily 0 group A member 3 Protein eagle Protein spready
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.