



Recombinant Enterobacteria phage HK022

Integrase (int)

Product Code	CSB-MP322283ECN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16407
Product Type	Recombinant Protein
Immunogen Species	Escherichia phage HK022 (Bacteriophage HK022)
Purity	≥85% (SDS-PAGE)
Sequence	MGRRRSHERR DLPPNLYIRN NGYYCYRDPR TGKEFGLGRD RRIAITEAIQ ANIELLSGNN RESLIDRIKG ADAITLHAWL DRYETILSER GIRPKTLLDY ASKIRAIRRK LPDKPLADIS TKEVAAMLNT YVAEGKSASA KLIRSTLV DV FREAI AEGHV ATNPVTATRT AKSEVRRSRL TANEYVAIYH AAEP LPIWLR LAMD LAVVTG QRVGDL CRMK WSDINDNHLH IEQSKTGAKL AIPLTLTIDA LNISLADTLQ QCREASSSET IIASKHHDPL SPKTVSKYFT KARNASGLSF DGNPPTFHEL RLSARLYRN QIGDKFAQRL LGHKSDSMAA RYRDSRGREW DKIEIDK
Source	Mammalian cell
Target Names	int
Protein Names	Recommended name: Integrase
Expression Region	1-357
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.