



Recombinant Human cytomegalovirus Protein UL24 (UL24)

Product Code	CSB-BP322311HWV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16760
Product Type	Recombinant Protein
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (HCMV)
Purity	≥85% (SDS-PAGE)
Sequence	YVRINAGQVL PVVWPPGWNL VLQEIETDED FKPEDVKAWS HYLCCQTRLA FVGRFVNEGV LSPDQQKKTA VCLISDEGYV FCYVREDTAV YYLARNLMEF ARVGLRAVET LHCMRYLTSS LVKRYFRPLL RAWSLGLDTM ARFIIRHHGQ FMPLTYPPGT ELRLCNLRCF ENSVEGGHLL RNIKTAFGMR VLGLGTVSLK GENAPFPHLR WPVDLIPIVV AYTGAVYACD VRDDRYIRVG DNLNTFMCLG LNLFFENRRF SGHNGIYDRV PDCPKGRQHR
Source	Baculovirus
Target Names	UL24
Protein Names	Recommended name: Protein UL24
Expression Region	21-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.