



Recombinant Burkholderia cepacia 2,2-dialkylglycine decarboxylase (dgdA)

Product Code	CSB-EP322324BXR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16932
Product Type	Recombinant Protein
Immunogen Species	Burkholderia cepacia (Pseudomonas cepacia)
Purity	>85% (SDS-PAGE)
Sequence	SLNDDATFW RNARQHLVRY GGTFFPMIIE RAKGSFVYDA DGRAILDFTS GQMSAVLGHC HPEIVSVIGE YAGKLDHLFS GMLSRPVVDL ATRLANITPP GLDRALLLST GAESNEAIR MAKLVTKYE IVGFAQSWHG MTGAAASATY SAGRKGVGPA AVGSFAIPAP FTYRPRFERN GAYDYLAELD YAFDLIDRQS SGNLAAFAIE PILSSGGIIE LPDGYMAALK RKCEARGMLL ILDEAQTGVG RTGTMFACQR DGVTPDILT SKTLGAGLPL AAVTSAIE ERAHELGYLF YTTHVSDPLP AAVGLRVLDV VQRDGLVARA NVMGDLRRG LLDLMERFDC IGDVRGRGLL LGVEIVKDRR TKEPADGLGA KITRECMNLG LSMNIVQLPG MGGVFRIAPP LTVSEDEIDL GLSLLGQAIE RAL
Source	E.coli
Target Names	dgdA
Protein Names	Recommended name: 2,2-dialkylglycine decarboxylase Short name= DGD EC= 4.1.1.64
Expression Region	2-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.