



Recombinant Human Zinc finger protein 30 (ZNF30)

Product Code	CSB-MP322332HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17039
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAHKYVGLQY HGSVTFEDVA IAFSQQEWES LDSSQRGLYR DVMLENYRNL VSMGHSRSPK HVALLEQWK EPEVTVRKDG RRWCTDLQLE DDTIGCKEMP TSENCPSFAL HQKISRQKPR ECQEYGKTLK QDSKPVQHER IHSSEKPNRC KECGKNFSNG HQLTIHQRLH VGEKPYKYEK CGKAFISGSA FVKHGRIHTG EKPLKCKQCG KTISGSYQLT VHKSIIHTGKK PYECGECGKA FLVYGKLTRH QSTHTGEKPF GCEECGKAFS TFSYLVQHQR IHTSEKPYEC KECGKAFSTS SPLAKHQRIH TGEKPYECKE CGKSFTVYGQ LTRHQSIHTG EKPFECKEKG KAFRLSSFLH AHQRIHAEIK PYGCKEKGRT FSRASYLVQH GRLHTGEKPY ECKEKGKAFS TGSYLVQHQR IHTGEKPYEC KECGKAFISR HQLTVHQRVH TGEKPYECKE CGKAFRVHVH LTQHRKIHTD VKPYECKEKG KTFSRASYLV QHSRIHTGKK PYECKEKGKA FSSGSYLVQH QRIHTGEKPY ECNKCGKAFT VYGQLIGHQS VHTGEKPFEC KECGKAFRLN SFLTEHQRVH TGEKPFKCKK CGKTFRYSSA LKVHLRKHMS VIP
Source	Mammalian cell
Target Names	ZNF30
Protein Names	Recommended name: Zinc finger protein 30 Alternative name(s): Zinc finger protein KOX28
Expression Region	1-623
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.