



# Recombinant Type-1 restriction enzyme Eco $\text{prrI}$ specificity protein (prrB)

<b>Product Code</b>	CSB-EP322344ENL
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P17222
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSELSYLEKL MDGVEVEWLP LSKVFNL RNG YTPSKTKKEF WANGDIPWFR MDDIRENGRI LGSSLQKISS CAVKGGKLP ENSILISTA TIGEHALITV PHLANQRFTC LALKESYADC FDIKFLFYC FSLAEWCRKN TTMSSFASVD MDGFKKFLIP RPCPDNPEKS LAIQSEIVRI LDKFSALTAE LTAELTAE LS MRKKQYNYR DQLLSFKEDE VEGKRKTLGE IMKMRAGQHI SAHNIIERKE ESYIYPCFGG NGIRGYVKEK SHDGEHLLIG RQGALCGNVQ RMKGQFYATE HAVVSVMPG INIDWAFHML TAMNLNQYAS KSAQPGLAVG KLQELKLFVP SIERQIYIAA ILDKFDLTN SITEVSRVKS SCARNSTNII EICYLVSRSR K
<b>Source</b>	E.coli
<b>Target Names</b>	prrB
<b>Protein Names</b>	Recommended name: Type-1 restriction enzyme Eco $\text{prrI}$ specificity protein Short name= S.Eco $\text{prrI}$ Alternative name(s): Type I restriction enzyme Eco $\text{prrI}$ specificity protein Short name= S protein
<b>Expression Region</b>	1-401
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.