



# Recombinant Type-1 restriction enzyme Ecoprrl specificity protein (prkB)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-YP322344ENL  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | P17222   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | MSELSYLEKL MDGVEVEWLP LSKVFNL RNG YTPSKTKKEF WANGDIPWFR<br>MDDIRENGRI LGSSLQKISS CAVKGGKLP ENSILISTA TIGEHALITV<br>PHLANQRFTC LALKESYADC FDIKFLFYC FSLAEWCRKN TTMSSFASVD<br>MDGFKKFLIP RPCPDNPEKS LAIQSEIVRI LDKFSALTAE LTAELTAE LS<br>MRKKQYNYR DQLLSFKEDE VEGKRKTLGE IMKMRAGQHI SAHNIIERKE<br>ESYIYPCFGG NGIRGYVKEK SHDGEHLLIG RQGALCGNVQ RMKGQFYATE<br>HAVVSVMPG INIDWAFHML TAMNLNQYAS KSAQPGLAVG KLQELKLFVP<br>SIERQIYIAA ILDKFDLTN SITEVSRVKS SCARNSTNII EICYLVSRSR K |
| <b>Source</b>            | Yeast  |
| <b>Target Names</b>      | prkB   |
| <b>Protein Names</b>     | Recommended name: Type-1 restriction enzyme Ecoprrl specificity protein<br>Short name= S.Ecoprrl Alternative name(s): Type I restriction enzyme Ecoprrl<br>specificity protein Short name= S protein   |
| <b>Expression Region</b> | 1-401  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |