



Recombinant *Saccharomyces cerevisiae* Protein KRE1 (KRE1)

| | |
|--------------------------|---|
| Product Code | CSB-YP322347SVG |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P17260 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast) |
| Purity | >85% (SDS-PAGE) |
| Sequence | AVTT QVTVVTVN VAG ALVTETTIWD PATAAAAATT TAQTGFFTTV FTTTNDVGTT VTLTQTVNRA TMLPTTTTST SSTGKTTTTV PTATSSLSSG LYLSTVTTTN DLGTTVTLTQ TFTHSSTSAT SSASSSVSSS VSSSGSSSSV KTTTSTGSAV AETGTRPDPS TDFTEPPVSA VTSLSIDSYI TITEGTTSTY TTTRAPSMW VTVVRQGNTI TVQTTFFVQRF SSQYVTVASP SVGSIGMGTL TGTVGVIKSA IKKTVSHN |
| Source | Yeast |
| Target Names | KRE1 |
| Protein Names | Recommended name: Protein KRE1 Alternative name(s): Killer toxin-resistance protein 1 |
| Expression Region | 27-288 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |