



# Recombinant Treponema pallidum Flagellar filament outer layer protein (flaA)

<b>Product Code</b>	CSB-YP322421TPR
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P18193
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Treponema pallidum (strain Nichols)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	DESVLIDFAK LNADIMADKS GGMTHNRRTV LDYASLADTS YTDEQKALMR SSLAVAQWEV VLNSSARNPV AHAASRVIEA PVSEGAKSFA GERVLGVRVL FPTWDSNANA MIKPAFVIPA YEVMAQVDDQ GNVQAPTEEE KASGKGRFED GYGVVKNVGV LKSIAVNTYG MNYPHGLYVM MRDQDGEVHR YFMGYLLFDS WKELVWNNPS YISDVSREV RLYPVYPAST PHVVFEGFMV TRDAAHAGGD YVG YFKDVKI IYDKAVLSTV RDFADEDLWG IQARREAERK RVEVARFGQQ QVLR YIEQEK LATEVGF TPS GGAQRQEEQQ
<b>Source</b>	Yeast
<b>Target Names</b>	flaA
<b>Protein Names</b>	Recommended name: Flagellar filament outer layer protein Alternative name(s): Sheath protein
<b>Expression Region</b>	21-350
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.