



Recombinant Escherichia coli Periplasmic protein trbC (trbC)

Product Code	CSB-MP322439ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P18473
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	SENVNTPEN RQFLKQQENL SRQLREKPDH QLKAWAEKQV LENPLQRSDN HFLDELVRKQ QASQDGKPRQ GALYFVSFSI PEEGLKRMLG ETRHFGIPAT LRGMVNNDLK TTAEAVLSLV KDGATDGVQI DPTLFSQYGI RTVPALVVFC SQGYDIIRGN LRVGQALEKV AATGDCRQVA HDLLAGKGDS GK
Source	Mammalian cell
Target Names	trbC
Protein Names	Recommended name: Periplasmic protein trbC
Expression Region	22-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.