



Recombinant Rat Cytohesin-1 (Cyth1)

Product Code	CSB-MP006467RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P97694
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MEDDDSYVPS DLTAEERQEL ENIRRRKQEL LADIQRLKEE IAEVANEIES LGSTEERKNM QRNKQVAMGR KKFNMDPKKG IQFLIENGLL KNTCEDIAQF LYKGEGLNKT AIGDYLGERD EFSIQVLHAF VELHEFTDLN LVQALRQFLW SFRLPGEAQK IDRMMEAFQA RYCQCNTGVF QSTDTCYVLS FAIIMLNTSL HNPNVKDKPT VERFIAMNRG INDGGDLPEE LLRNLYESIK NEPFKIPEDD GNDLTHFFN PDREGWLLKL GGGRVKTWKR RWFILTDNCL YFEYTTDKE PRGIIPLENL SIREVEDSKK PNCFELYIPD NKDQVIKACK TEADGRVVEG NHTVYRISAP TPEEKEDWIK CIKAAISRDP FYEMLAARKK KVSSTKRH
Source	Mammalian cell
Target Names	Cyth1
Protein Names	Recommended name: Cytohesin-1 Alternative name(s): PH, SEC7 and coiled-coil domain-containing protein 1 SEC7 homolog A Short name= rSec7-1
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein is a member of the PSCD family. Members of this family have identical structural organization that consists of an N-terminal coiled-coil motif, a central Sec7 domain, and a C-terminal pleckstrin homology (PH) domain. The coiled-coil motif is involved in homodimerization, the Sec7 domain contains guanine-nucleotide exchange protein (GEP) activity, and the PH domain interacts with phospholipids and is responsible for association of PSCDs with membranes. Members of this family appear to mediate the regulation of protein sorting and membrane trafficking. This gene is highly expressed in natural killer and peripheral T cells, and regulates the adhesiveness of integrins at the plasma membrane of lymphocytes. The encoded protein is 83% homologous to that of CYTH2.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.