



# Recombinant Enterobacteria phage T4 Baseplate structural protein Gp8 (8)

<b>Product Code</b>	CSB-BP322492EDZ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P19062
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Enterobacteria phage T4 (Bacteriophage T4)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNDSSVIYRA IVTSKFRTEK MLNFYNSIGS GPKNTIFIT FGRSEPWSSN ENEVGFAPPY PTDSVLGVTD MWTMMGTVK VLPSMLDAVI PRRDWGDTRY PDPYTFRIND IVVCNSAPYN ATEGAGWLV YRCLDVPDTG MCSIASLTDK DECLKLGKWK TPSARSMTPP EGRGDAEGTI EPGDGYVWEY LFEIPPDVSI NRCTNEYIVV PWPEELKEDP TRWGYEDNLT WQQDDFGLIY RVKANTIRFK AYLDSVYFPE AALPGNKGFR QISIITNPLE AKAHPNDPNV KAEKDYYDPE DLMRHS GEMI YMENRPPIIM AMDQTEEINI LFTF
<b>Source</b>	Baculovirus
<b>Target Names</b>	8
<b>Protein Names</b>	Recommended name: Baseplate structural protein Gp8 Alternative name(s): Baseplate wedge protein 8
<b>Expression Region</b>	1-334
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.