



Recombinant Enterobacteria phage T4 Baseplate structural protein Gp8 (8)

Product Code	CSB-MP322492EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19062
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MNDSSVIYRA IVTSKFRTEK MLNFYNSIGS GPKNTIFIT FGRSEPWSSN ENEVGFAPPY PTDSVLGVTD MWTHMMGTVK VLPSMLDAVI PRRDWGDTRY PDPYTFRIND IVVCNSAPYN ATEGAGWLV YRCLDVPDTG MCSIASLTDK DECLKLGKWK TPSARSMTPP EGRGDAEGTI EPGDGYVWEY LFEIPPDVSI NRCTNEYIVV PWPEELKEDP TRWGYEDNLT WQQDDFGLIY RVKANTIRFK AYLDSVYFPE AALPGNKGFR QISIITNPLE AKAHPNDPNV KAEKDYYDPE DLMRHS GEMI YMENRPPIIM AMDQTEEINI LFTF
Source	Mammalian cell
Target Names	8
Protein Names	Recommended name: Baseplate structural protein Gp8 Alternative name(s): Baseplate wedge protein 8
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.