



Recombinant Simian immunodeficiency virus Protein Nef (nef)

| | |
|--------------------------|---|
| Product Code | CSB-MP322539SHD |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P19501 |
| Product Type | Recombinant Protein |
| Immunogen Species | Simian immunodeficiency virus (isolate PBJ14/BCL-3) (SIV-sm) (Simian immunodeficiency virus sooty mangabey monkey) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | GGVTSKKQR RRGGLYERL LQARGETYGR LWEGLEGEYS QSQDASGKGL SSLSCEPQKY CEGQFMNTPW RNPATERAKL DYRQQNMDDV DSADLVGCPV SPRVPVRIMT YKLAIMSHF IKEKGGLEGI YYSDRRHKIL DLYLEKEEGI IPDWQNYTAG PGRIYPMFFG WLWKLVPVNV SDEAQEETH YLMHPAQTSQ WDDPWGEVLA WKFDPKLAYN YKAFVEHPEE FGSQSGLSKE EVQRRLTARG LLKMADKKKT S |
| Source | Mammalian cell |
| Target Names | nef |
| Protein Names | Recommended name: Protein Nef Alternative name(s): 3'ORF Negative factor Short name= F-protein |
| Expression Region | 2-261 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |