



# Recombinant Type-2 restriction enzyme BanI (banIR)

<b>Product Code</b>	CSB-MP322573BAF
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P19887
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Aneurinibacillus aneurinilyticus (Bacillus aneurinolyticus)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AQLKYNKDI DELERNAAKW WPDFLAKKES STSIIPKLVE SQDAFISLLN LSKNNPFDIF QLIDASKFPP NLFLKHLVVL TDFGGEPLNR LNQNFDLSLFP MIPYGIHYIT KVLGKFEFFW NEKKYEYVFQ ELPVTSLTNS KLKIDGASIS KTVPLSDLYK DVIVLLMFGA NAVNSEVSEV LMKCEVGNLI GKTDELKKFI KERYIFVSRI TGGAEANTLG QVAQTHVIDF LRTRFGSKGH DIKSNGHIEG VTHNDGQTLT TFDVVIKKG SVAIEISFQ VTTNSTIERK AGQAKARYDM VSDTGNYIAY IIDGAGNFQR KNAITTCNN SHCTVAYTEE ELNVLLKFIL EKLE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	banIR
<b>Protein Names</b>	Recommended name: Type-2 restriction enzyme BanI Short name= R.BanI EC= 3.1.21.4 Alternative name(s): Endonuclease BanI Type II restriction enzyme BanI
<b>Expression Region</b>	2-354
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.