



Recombinant Neisseria gonorrhoeae Major outer membrane protein P.IB(porB)

Product Code	CSB-BP322601NFZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20148
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Neisseria gonorrhoeae
Purity	>85% (SDS-PAGE)
Sequence	D VTLYGAIKAG VQTYRSVEHT DGKVSQVETG SEIADFGSKI GFKGQEDLGN GLKAVWQLEQ GASVAGTNTG WGNKQSFVGL KGGFGTIRAG SLNSPLKNTG ANVNAWESGK FTGNVLEISG MAQREHRYLS VRYDSPEFAG FSGSVQYAPK DNSGSNGESY HVGLNYQNSG FFAQYAGLFQ RYGEGTKKIE YDDQTY SIPS LFVEKLQVHR LVGGYDNNAL YVSVAQQQD AKLYGAMSGN SHNSQTEVAA TAAYRFGNVT PRVSYAHGFK GTVDSANHDN TYDQVVVGAE YDFSKRTSAL VSAGWLQGGK GADKIVSTAS AVVLRHKF
Source	Baculovirus
Target Names	porB
Protein Names	Recommended name: Major outer membrane protein P.IB Short name= PIB Short name= Protein IBA Alternative name(s): Porin
Expression Region	20-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.