



Recombinant Solanum tuberosum Cysteine protease inhibitor 1

Product Code	CSB-EP322630FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20347
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	LVLPEVYD QDGNPLRIGE RYIINNPLLG AGAVYLYNIG NLQCPNAVLRQ HMSIPQFLGE GTPVVFVRKS ESDYGDVVRV MTVVYIKFFV KTTKLCVDQT VWKVNDEQLV VTGGKVGREN DIFKIMKTDL VTPGGSKYVY KLLHCPSHLG CKNIGGNFKN GYPRLVTVDD DKDFIPFVFI KA
Source	E.coli
Protein Names	Recommended name: Cysteine protease inhibitor 1 Alternative name(s): P340 P34021 PCPI 8.3
Expression Region	43-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.