




Recombinant Moraxella bovis Fimbrial protein 1

Product Code	CSB-EP322675MHZ
Abbreviation	Recombinant Moraxella bovis Fimbrial protein 1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20657
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Moraxella bovis
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	FTLIELMIVIAIIGILAAIALPAYQDYISKSQTTRVSGELAAGKTAVDAAALFEGKTPV LSEESSTSKENIGLTSSETSTKPRSNLMASVELTGFADNGAGTISATLGNKANK DIAKTIVITQERTTDGVWTCKIDGSQAQYKEKFNPTGCVKK
Research Area	Others
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Fimbrial protein 1Alternative name(s): Alpha-pilin Fimbrial protein I I pilin
Expression Region	7-159aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	29.0 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.