



# Recombinant *Oryza sativa* subsp. japonica Glucose-1-phosphate adenylyltransferase small subunit, chloroplastic/amyloplastic (AGPS)

<b>Product Code</b>	CSB-YP322766OFG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P15280
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AVSDSK SSQTCLDPDA STSVLGIILG GGAGTRLYPL TKKRAKPAVP LGANYRLIDI PVSNCLNSNI SKIYVLTQFN SASLNRHLSR AYGNNIGGYK NEGFVEVLAA QQSPDNPWF QGTADAVRQY LWLFEEHNVM EFLILAGDHL YRMDYEKFIQ AHRETSDIT VAALPMDEKR ATAFGLMKID EEGRIVEFAE KPKGEQLKAM MVDTTILGLD DVRAKEMPYI ASMGIYVISK NVMLQLLREQ FPGANDFGSE VIPGATNIGM RVQAYLYDGY WEDIGTIEAF YNANLGITKK PVPDFSFYDR SAPIYTQPRH LPPSKVLDAD VTDSVIGEGC VIKNCKIHHS VVGLRSCISE GAIIEDSLM GADYYETead KKLLGEKGGI PIGIGKNCHI RRAIDKNAR IGDNVKIINV DNVQEAARET DGYFIKSGIV TVIKDALLPS GTVI
<b>Source</b>	Yeast
<b>Target Names</b>	AGPS2
<b>Protein Names</b>	Recommended name: Glucose-1-phosphate adenylyltransferase small subunit, chloroplastic/amyloplastic EC= 2.7.7.27 Alternative name(s): ADP-glucose pyrophosphorylase ADP-glucose synthase AGPase B Alpha-D-glucose-1-phosphate adenylyl
<b>Expression Region</b>	65-514
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.