



Recombinant Drosophila melanogaster Protein henna (Hn)

Product Code	CSB-MP322921DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17276
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MYQRQVSFDK PTRVEDSAYI VEGVDIKEAR NTCLLFSPKD SSLSSGALAN ILAIFKKHDI NLVHIESRSS LRVPGYEFFV EADGKSGALG KAIEDVKEQC SYFNIISRDI KDNATAVPWF PRRIRDLDRF ANQILSYGSE LDADHPGFTD PEYRKRRKYF ADIAYNYKHG EPLPHVDYTK EEIETWGIIF RNLTCLYKTH ACREYNHVFP LLVDNCGFRE DNIPQLEDVS NFLRDCTGFT LRPVAGLLSS RDFLAGLAFR VFHSTQYIRH PSKPMYTPEP DVCHELMGHV PLFADPAFAQ FSQEIGLASL GAPDDYIEKL STIFWFTVEY GVCRRQEGELK AYGAGLLSSY GELEYCLTDK PQLKDFEPEV TGVTKYPITQ FQPLYVADS FETAKEKTIK FANSIPRPFV VRYNAYTQSV EVLDSKPQIS NLMDNINSEF QILQNAVAKL RV
Source	Mammalian cell
Target Names	Hn
Protein Names	Recommended name: Protein henna EC= 1.14.16.1 EC= 1.14.16.4 Alternative name(s): Phe-4-monooxygenase Phenylalanine-4-hydroxylase Short name= PAH Tryptophan 5-hydroxylase Short name= TRH Tryptophan 5-monooxyg
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.