



# Recombinant Lactobacillus paracasei D-2-hydroxyacid dehydrogenase (NAD<sup>+</sup>)

<b>Product Code</b>	CSB-YP322950LAY
<b>Abbreviation</b>	Recombinant Lactobacillus paracasei D-2-hydroxyacid dehydrogenase (NAD <sup>+</sup> )
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P17584
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Lactobacillus paracasei
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKIIAYGARV DEIQYFKQWA KDTGNTLEYH TEFLDENTVE WAKGFDGINS LQTPPYAAGV FEKMHAYGIK FLTIRNVGTD NIDMTAMKQY GIRLSNPAY SPAAIAEFAL TDTLYLLRNM GKVQAQLQAG DYEKAGTFIG KELGQQTGVG MGTGHIGQVA IKLFKGFQAK VIAYDPYPMK GDHPDFDYVS LEDLFKQSDV IDLHVPGIEQ NTHIINEAAF NLMKPGAIVI NTARPNLIDT QAMLSNLKSG KLAGVGIDTY EYETEDLLNL AKHGSFKDPL WDELLGMPNV VLSPHIAYYT ETAVHNMVYF SLQHLVDFLT KGETSTEVTG PAK
<b>Source</b>	Yeast
<b>Target Names</b>	N/A
<b>Protein Names</b>	Recommended name: D-2-hydroxyisocaproate dehydrogenase Short name= D-HICDH EC= 1.1.1.-
<b>Expression Region</b>	1-333
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.