



Recombinant Enterobacteria phage T4 Protein inh (inh)

Product Code	CSB-MP322987EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P18058
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MIDKDYIAEL KALDDNKEAK AKLAEYAEQF GIKVKKNKSF DNIVVDIEEA LQKLASEPMP ETDGLSIKDL IDAADA A EGL KYDDEEVNPE AALLIDSPIK SDIKIEVVET DKIPENTDVL IEDTPFVEEK FEQAVAEIIE SEKPSVFTLP ENFSPNLQLI GKNLGFCTVP WWIYQWIAET PDWKSHPTSF EHASAHQTLF SLIYYINRDG SVLIRETRNS SFVTLK
Source	Mammalian cell
Target Names	inh
Protein Names	Recommended name: Protein inh Alternative name(s): Inhibitor protein
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.