



Recombinant Cytolysin (vvhA)

Product Code	CSB-MP323085VFI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19247
Product Type	Recombinant Protein
Immunogen Species	Vibrio vulnificus (strain CMCP6)
Purity	≥85% (SDS-PAGE)
Sequence	QEYVPIVEKP IYITSSKIKC VLHTSGDFNA TRDWCNAGAS IDVRVNVAQM RSVQSATSDG FTPDAKIVRF TVDADKPGTG IHLVNELQQD HSWFQSWANR RTYIGPFASS YDLWVKPVSG YTPKKARDLP QNENKNYQHR DTYGYSIGIN GKVGAEVNKD GPKVGGEVSG SFTYNYSKTL VFDTKDYRIN NRSSLSDFDI SFEREFGECD ELRRQELGCY FTAAHWGSGW VFDKTKFNPI SYSNFKPNYD VLYEAPVSET GVTDFEMGVK LNYRARFGTV LPSALFSVYG SAGSSTNSST VKQRIRIDWN HPLFEAEAHV TLQSLSNNDL CLDVYGENGD KTVAGGSVNG WSCHGSWNQV WGLDKEERYR SRVASDRCLT VNADKTLTVE QCGANLAQKW YWEGDKLISR YVDGNNTRYL LNIVGGRNVQ VTPENEANQA RWKPTLQQVK L
Source	Mammalian cell
Target Names	vvhA
Protein Names	Recommended name: Cytolysin
Expression Region	21-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.