



Recombinant Salmonella typhimurium Maltose-binding periplasmic protein (malE)

Product Code	CSB-YP323126SXB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19576
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	KIEE GKLVIWINGD KGYNGLAEVG KKFEQDTGIK VTVEHPDKLE EKFPQVAATG DGPDIIFWAH DRFGGYAQSG LLAEVTPDKA FQDKLYPFTW DAVRYNGKLI AYPIAVEALS LIYNKDLVPN PPKTWEEIPA LDKELKVKGK SAIMFNLQEP YFTWPLIAAD GGYAFKFENG KYDVKDVGVN NAGAKAGLTF LIDMIKNKNM SADTDYSIAE AAFNKGETAM TINGPWAWSN IDKSKVNYGV TLLPTFKGKP SKPFVGVLSA GINAASPNKE LAKEFLENYL LTDQGLEAVN KDKPLGAVAL KSFQEQLAKD PRIAATMDNA QKGEIMPNIQ QMSAFWYAVR TAVINAASGR QTVDAALKDA QSRITK
Source	Yeast
Target Names	malE
Protein Names	Recommended name: Maltose-binding periplasmic protein Alternative name(s): MBP MMBP Maltodextrin-binding protein
Expression Region	27-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.