



Recombinant Protein lin-11 (lin-11)

Product Code	CSB-EP323195CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20154
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MHSSSSFIIT SLEEEEEKPP AHHLHQQSIE DVGSVTSSAT LLLLDSATWM MPSSTTQPHI SEISGNECAA CAQPILDYV FTVLGKCVHW SCLRCCDCRA PMSMTCFSRD GLILCKTDFS RRYSQRCAGC DGKLEKEDLV RRARDKVFHI RCFQCSVCQR LLDTGDQLYI MEGNRFVCQS DFQTATKTST PTSIHRPVSN GSECNSDVEE DNVDACDEVG LDDGEGDCGK DNSDDNSAK RRGPRTTIKA KQLETLKNAF AATPKPTRHI REQLAAETGL NMRVIQVWFQ NRRSKERRMK QLRFGGYRQS RRPRRDDIVD MFPNDQQFYF PPPPSNVQFF CDPYTTSPNN PETIQMAPQF AVPTENMMNV PEPYTEQSAT PPEFNEDTFA CIYSTDLGKP TPVSW
Source	E.coli
Target Names	lin-11
Protein Names	Recommended name: Protein lin-11 Alternative name(s): Abnormal cell lineage protein 11
Expression Region	1-405
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.